

Medical Profile Setup

Cycle 1

11/09/2020

Test 1

Screen Dimension

- 6.1-inch (diagonal)
- 5.94 inch (height)
- 2.98 inches (width)

Model: iPhone 11

id	Description	Steps to reproduce	Expected Result	Actual Result	Date	Name	Result
1	Verify that all input fields are shown and fitted to screen, along with the 'Show in Lock Screen' CheckBox and 'Save' button	Navigate to Medical Profile Setup Screen	Medical Form is displayed with all fields fitted to screen, along with the bottom checkbox and 'Save' button.	Medical Form is displayed with all fields, along with the bottom checkbox and 'Save' button. But is NOT fitted to screen, showing a bottom white margin.	11/04/20	Hannah	FAIL
2	Verify that Keyboard does not cover text input boxes upon filling out.	Navigate to Medical Profile Setup Screen	Medical Form is displayed with all fields fitted to screen, along with the bottom checkbox and 'Save' button.	Medical Form is displayed with all fields, along with the bottom checkbox and 'Save' button. But is NOT fitted to screen, showing a bottom white margin.	11/04/20	Hannah	FAIL

		Click on any of the input text fields	Keyboard is pushed under that current text field.	Keyboard is NOT pushed under that current text field, ultimately hiding the current box completely.			
--	--	---------------------------------------	---	---	--	--	--

<div><div>Id: 1</div><div><div><div>4:06</div><div>✕</div></div><div><div>Medical Information</div><div><div>Age</div><div>Sex</div><div>Doctors/Caregiver Name</div><div>Next of kin</div><div>Blood Type</div><div>Medical Conditions</div><div>Allergies</div><div>Medication List</div><div>Additional Information</div><div><div><input type="radio"/> Show in locked screen</div><div>Save</div></div></div></div></div></div>	<div><div>Id: 2</div><div><div><div>17:07</div><div>✕</div></div><div><div>Medical Information</div><div><div>Age</div><div>Sex</div><div>Doctors/Caregiver Name</div><div>Next of kin</div><div>Blood Type</div><div>Additional Information</div></div></div><div><div>Hahahahaha I Hey</div><div><div>QWERTYUIOP</div><div>ASDFGHJKL</div><div>↑ZXCVBNM↵</div><div>123spacereturn</div><div>😊🎤</div></div></div></div></div>
--	--

Test 2

Screen Dimension

- 4.7 inch (diagonal)
- 5.45 inches inch (height)
- 2.65 inches (width)

Model: iPhone 8

id	Description	Steps to reproduce	Expected Result	Actual Result	Date	Name	Result
1	Verify that all input fields are shown and fitted to screen, along with the 'Show in Lock Screen' CheckBox and 'Save' button	Navigate to Medical Profile Setup Screen	Medical Form is displayed with all fields fitted to screen, along with the bottom checkbox and 'Save' button.	Medical Form is displayed with all fields, along with the bottom checkbox and 'Save' button. But it is NOT fitted to the screen, the user needs to scroll the page.	11/04/20	Hannah	FAIL
2	Verify that Keyboard does not cover text input boxes upon filling out.	Navigate to Medical Profile Setup Screen	Medical Form is displayed with all fields fitted to screen, along with the bottom checkbox and 'Save' button.	Medical Form is displayed with all fields, along with the bottom checkbox and 'Save' button. But it is NOT fitted to the screen, the user needs to scroll the page.	11/04/20	Hannah	FAIL
		Click on any of the input text fields	Keyboard is pushed under that current text field.	Keyboard is NOT pushed under that current text field, ultimately hiding the current box			

				completely.			
--	--	--	--	-------------	--	--	--

Id: 1

CarrierMedical Information4:20 PM

Age

Sex

Doctors/Caregiver Name

Next of kin

Blood Type

Medical Conditions

Allergies

Medication List

Additional Information

Show in locked screen

Save

Id: 2

CarrierMedical Information4:20 PM

Age

Sex

Doctors/Caregiver Name

Next of kin

Blood Type

Medical Conditions

ITheI'm

QWERTYUIOP

ASDFGHJKL

↑ZXCVBNM↵

123😊🎤space↵

Test 3

Screen Dimension

- 4.0 inch (diagonal)
- 4.87 inches inch (height)
- 2.31 inches (width)

Model: iPhone 5

id	Description	Steps to reproduce	Expected Result	Actual Result	Date	Name	Result
1	Verify that all input fields are shown and fitted to screen, along with the 'Show in Lock Screen' CheckBox and 'Save' button	Navigate to Medical Profile Setup Screen	Medical Form is displayed with all fields fitted to screen, along with the bottom checkbox and 'Save' button.	Medical Form is displayed with all fields, along with the bottom checkbox and 'Save' button. But it is NOT fitted to the screen, the user needs to scroll the page.	11/04/20	Hannah	FAIL
2	Verify that Keyboard does not cover text input boxes upon filling out.	Navigate to Medical Profile Setup Screen	Medical Form is displayed with all fields fitted to screen, along with the bottom checkbox and 'Save' button.	Medical Form is displayed with all fields, along with the bottom checkbox and 'Save' button. But it is NOT fitted to the screen, the user needs to scroll the page.	11/04/20	Hannah	FAIL
		Click on any of the input text fields	Keyboard is pushed under that current text field.	Keyboard is NOT pushed under that current text field, ultimately hiding the current box completely.			

Id: 1	Id: 2
<div data-bbox="215 315 655 1099"><div>Carrier 4:12 PM</div><div>×</div><div>Medical Information</div><div>Age</div><div>Sex</div><div>Doctors/Caregiver Name</div><div>Next of kin</div><div>Blood Type</div><div>Medical Conditions</div><div>Allergies</div></div> <div data-bbox="215 1111 655 1899"><div>Carrier 4:12 PM</div><div>Sex</div><div>Doctors/Caregiver Name</div><div>Next of kin</div><div>Blood Type</div><div>Medical Conditions</div><div>Allergies</div><div>Medication List</div><div>Additional Information</div><div><input type="radio"/> Show in locked screen</div><div>Save</div></div>	<div data-bbox="817 315 1385 831"><div>Carrier 5:23 PM</div><div>Sex</div><div>Doctors/Caregiver Name</div><div>Next of kin</div><div>Blood Type</div><div>Medical Conditions</div></div> <div data-bbox="817 875 1385 1328"><div>I The I'm</div><div>Q W E R T Y U I O P</div><div>A S D F G H J K L</div><div>↑ Z X C V B N M ↵</div><div>123 😊 space return</div></div>